The High Level Disinfectant of Choice in...

- ICU
- SICU
- ICCU
- NICU

- Cath Lab
- Dialysis Unit
- Dental Clinic
- Operation Theatre

- Endoscopy
- **Blood Bank**
- Microbiology Lab
- Emergency Ward

USAGE

- ▶ Critical Areas: (1:64) 150 ml in 10 Litre of water for mopping and wiping walls, floors, furnitures, operation tables, O.T. lamps & other equipments.
- Non-Critical Areas: (1:128) 75 ml in 10 Litre of water for mopping in corridors, waiting areas, reception areas & hospital wards.
- ▶ Fogging: 15 ml of GnD -125 + concentrate with 1 Litre tap water for every 1000 cubic feet in critical or high risk areas like ICU, ICCU, NICU, Cath Lab, AKD, O.T., Emergency wards etc.

COMPARATIVE EFFICIENCY					
	GnD-125+	Brands (Containing Formaldehyde/Hydrogen peroxide & Silver nitrate, Dodecyl benzene-sulfonate, Sulfamic acid etc.)			
Contact time	10 min.	≥10 min.			
 Dilution for fogging 	15 ml	≥100 ml			
 Post fogging- Re-entry time 	45 min.	≥1 hr.			
 Activity in hard water 	Upto 800 ppm	X			
Soil tolerance	95%	5%			
Eco-friendly		X			
Non-irritant		X			
Non-corrosive	V	X			
 Non-carcinogenic 		X			
Non-staining		x			
Non-sticky		X			
Time consuming	X				

100% Sterilization! 100% Secure! 100% Safe!

ACCREDITATIONS

















Nanz Med Science Pharma Pvt. Ltd.

(A Subsidiary of Med Science Canada Inc.)





Third Generation Twin Chain QAC Disinfectant for Surface & Environment



Super fast & Super active Disinfectant than ever

GnD-125+

- N-Alkyl Dimethyl Benzyl Ammonium Chloride (C12-C18) 2.37%w/v
- N-Alkyl Dimethyl Ethyl Benzyl Ammonium Chloride (C12-C14) 2.37%w/v





- Eradicates more than 125 types of pathogens like H1N1, H5N1, HIV, HCV, HBV, Corona Virus, Polio etc.
- No cross resistance reported
- Non-corrosive & Non-irritating
- Non-carcinogenic & Non-hazardous
- Ideal for fogging to disinfect critical
 & non-critical areas
- Highly biodegradable

Packing:



200 ml

SHELF LIFE

3 YEARS (PACKED) 64 DAYS (DILUTED SOLUTION)

Description

- GnD-125 + is a third generation twin chain Quaternary Ammonium Compounds
- 99.99% Accuracy against 125 + types of pathogens including H1N1, H5N1, HBV, HCV, HIV, Corona Virus, Polio etc.
- A highly stabilized twin chain QAC disinfectant solution for hard surfaces
- Ideal for various types of medical & healthcare facilities
- Absolutely no report of resistance buildup

Results of plate exposure method showing the number of bacterial & fungal colonies before and after sterilisation with third generation quaternaries.

Int J Med Res Rev 2015;3(3):273-277

Operation Theatre	Number Bacterial Colonies		
	Before Sterilisation	After Sterilisation	
General Surgery	4	0	
ENT	3	0	
OBG	5	0	
Ophthalmology	4	0	
Orthopaedics	2	0	

Operation Theatre	Number of Fungal Colonies		
	Before Sterilisation	After Sterilisation	
General Surgery	8	0	
ENT	3	0	
OBG	6	0	
Ophthalmology	9	0	
Orthopaedics	18	0	

The study has shown that the quaternary ammonium compounds are effective disinfectant and sterilisation agent for operation theatres.

Int J Med Res Rev2015;3(3):273-277

The results showed the complete removal of micro organism after 45 minutes of fogging by third generation Quaternaries.

Int J Med Res Rev2015;3(3):273-277

The Potent **Disinfectant**

which is highly Environment Friendly!!