

# Germicidal Liquid Cleaner For Industrial Equipments

## Nanzilon-HC

Chlorhexidine Gluconate Solution 7.5% v/v,  
Cetrimide 15% w/v, Iso propyl Alcohol 6-8% v/v

Hospital  
Concentrate



1 Litre



500 ml

### Salient Features:

- Strong germicidal & detergent property
- For cleaning of Industrial equipments
- Swabbing floors, utility & side tables
- Instant cleaning with non-staining residue

### How to Use:

- **For industrial linen and furnitures:**  
1 in 100 part aqueous (100 ml diluted to 1 Litre water)
- **For Instruments:**  
10ml diluted to 1 Litre water & add sodium nitrate (4gm/Litre) to prevent rusting
- **For rubber appliances:**  
1 in 30 parts of aqueous (35 ml diluted to 1 Litre water)

### Common Usage Areas:

- In laundries, In working areas for cleaning equipments, Swabbing floors, Utility & side tables



Also available in pack sizes

\*20 Litre Pack is also available